

Developing a high-quality, sample-to-result Hereditary Breast and Ovarian Cancer panel assay pipeline for a novel sequencing platform

A. Bhattacharjee¹, G. Liu¹, E. Frise¹, S. Kruglyak², B. Lajoie², K. Blease², J. Zhao², M. Reese¹, S. Levy².
¹Fabric Genomics, Oakland, CA; ²Element Biosciences, San Diego, CA

INTRODUCTION

Hereditary Breast and Ovarian Cancer (HBOC) testing has remained inadequate, plagued by medical history-based testing, high test cost, and challenges in NGS data quality. This has led to a model of centralized testing on high-throughput sequencers that always do not fit well with independent clinical testing. We explore herein alternative solutions and improvement approaches.

METHODS

- Compared 5 Coriell cell line and 46 patient samples across two short read sequencing platforms (Element and Illumina NextSeq550)
- Used TruSight Hereditary Cancer Panel (113 genes) that contain HBOC Genes to compare performance on Illumina platform for SNV and CNV calls to identify improvement areas

RESULTS

- Target read depth highly correlated across the two sequencer platforms
- High concordance across clinically significant variants: 277/279 sites
- Benchmark Study: Performance Metrics on Illumina Sequencer: High Sensitivity [TP/(TP+FP)] and PPV [TP/(TP+FP)]/Specificity [TN/(TN+FP)]
- Improved FG Panel design for CNV detection

20,406 total variants in 51 samples with high concordance between platforms

Clinically Classified Variants by FG	Concordant Genotypes (Element/Illmn)	Discordant Genotypes (Element/Illmn)
279	277	2

Only two discordant variants. Both located at sites next to a long homopolymer, where read mapping can be difficult and classified as benign variants

Performance on Illumina Sequencer

Variant Type	Reference Material	Performance
SNVs and small Indels in Genome-in-a-Bottle consortium +/- 10-bp	NIST ID HG001, HG002, HG003, HG004 (Coriell ID NA12878, NA24385, NA24149, NA24143), across 20 HBOC genes (see above)	<ul style="list-style-type: none">Sensitivity: 99.38% for HBOC and 98.49% for 113 gene panelPPV HBOC or panel: ~99%
SNVs and small Indels, Challenging variants	Seracare Seraseq™ Inherited Cancer DNA Mix v1 of 23 variants across (BRCA1, BRCA2, CHEK2, EPCAM, MLH1, MSH2, MSH6, PMS2)	Sensitivity: 100% with the exception of 2 challenging indel variants-126-bp insertion and a 300-bp AluY insertion in BRCA2 that require visual review
Large and Small CNVs	Orthogonally validated collection of cases from (a) 33 Coriell samples spanning large and small CNVs. (b) ICR 96 dataset for 56 CNVs across 14 genes (ATM, BRCA1, BRCA2, CHEK2, EPCAM, MLH1, MSH2, MSH6, NF1, PALB2, PMS2, PTEN, RAD51C, TP53)	<ul style="list-style-type: none">Sensitivity: ~97% for HBOC; ~95% for 113 gene panelSpecificity: >99.98 %


REFERENCES

1. Mahamdallie S, Ruark E, Yost S et al. The ICR96 exon CNV validation series: a resource for orthogonal assessment of exon CNV calling in NGS data. Wellcome Open Res 2017, 2:35 (doi.org/10.12688/wellcomeopenres.11689.1)

Samples used in the study

	Sample ID	Description	Size	Gene
5 Coriell	1 NA14091	5382insC (exon20)	1 bp	BRCA1
	2 NA11630	317ins5	5 bp	MEN1
	3 NA11410	3149delC (resulting in truncation of gene product)	1 bp	APC
	4 NA10080	781C>T (Gln261Ter)	1 bp	PTEN
	5 NA16533	19-bp deletion in exon 2	19 bp	CDKN2A
46 patients 51 samples	6 Patient 1			
	7 Patient 2			
	...			
	51 Patient 46			


Sequencing Platform Workflows



BRCA1 mapped reads (IGV)



CNV Calling on Illumina TruSight/Illumina Sequencer




FG PAN CANCER panel (113 genes)

HBOC PANEL		
Breast Cancer	Ovarian Cancer	Breast & Ovarian Cancer (HBOC)
ATM	ATM	ATM
BARD1	BRCA1	BARD1
BRCA1	BRCA1	BRCA1
BRCA2	BRCA2	BRCA2
BRIP1	BRIP1	BRIP1
CHEK2	CHEK2	CHEK2
CDH1	EPCAM*	EPCAM*
MLH1	MLH1	MLH1
MSH2	MSH2	MSH2
MSH6	MSH6	MSH6
NF1	NF1	NF1
PALB2	PALB2	PALB2
PMS2*	PMS2*	PMS2*
PTEN	PTEN	PTEN
RAD51C	RAD51C	RAD51C
RAD51D	RAD51D	RAD51D
STK11	STK11	STK11
TP53	TP53	TP53

■ NCCN v1.2022 Risk and Management, Evidence
* NCCN v1.2022 Risk and Management, Limited Evidence

- Colorectal
- Pancreatic
- Melanoma
- Pediatric
- Brain

Correlation: Read depth per target




DISCUSSION

- Improvements in CNV calling (secondary) and panel design is expected to boost performance (high sensitivity) in short read sequencing platforms

CONCLUSION

- Combination of joint CNV calling and panel design is expected to provide a high-quality short-read sequencing solution for HBOC and pan cancer testing

FG Apollo (CNV exon) Performance



Panel Probe Design for FG Panel vs Illumina TruSight (BRCA1)

